


FD4000

TECHNICAL SPECIFICATIONS



SENSOR PLATE

QUANTITY	2 per system
DIMENSIONS	350 mm W x 700 mm L x 78 mm H
WEIGHT	23 kg (46 kg per pair)
SAMPLE RATE	up to 1000 Hz
CAPACITY	2,000 kg
SAFE OVERLOAD	3,000 kg
MAXIMUM OVERLOAD	4,000 kg
MEASURING RANGE	0 - 1,000 kg (Fz only)
RESOLUTION	c. 15 g / 0.15 N
LINEARITY AND HYSTERESIS	<0.05% over calibrated range

ALUMINIUM FRAME

DIMENSIONS	820 mm W x 820 mm L x 78 mm H
WEIGHT	10 kg

ELECTRICAL

POWER	Output 24V from supplied 90-264V PSU
CONNECTIVITY	USB

SOFTWARE SYSTEM REQUIREMENTS

PROCESSOR	Intel i5 / i7
OPERATING SYSTEM	Windows x64 version 7 or later
RAM	8GB
HARD DRIVE	Minimum 500MB for install
PORTS	1 x USB 2.0 / USB 3.0